

## CAD4TB Performance Summary

Delft Imaging's AI-powered CAD4TB is the most rigorously validated CAD technology for TB screening and triaging. Over 100 publications highlight its technicality, operational performance, cost-effectiveness and reliability.

Various studies suggest that CAD4TB performs significantly better than human readers and meets the WHO's Target Product Profile, reduces Xpert MTB/RIF cartridges consumption, and increases daily screening throughput, thus contributing to reducing cost and making TB services more accessible, accurate, efficient and effective.

### WHO Target Product Profile

CAD4TB is one of the three CAD software packages included in the evaluation to support WHO's recommendation on using CAD for TB screening and triage for individuals aged 15 years and older<sup>1</sup> (meeting the TPP accuracy standard of 90% sensitivity and 70% specificity).

An independent, head-to-head comparison of five CAD algorithms accuracy for presumptive TB, including key groups (PLWH and diabetes), across seven countries (India, Madagascar, South Africa, Tanzania, the Philippines, South Africa, and Vietnam) concluded that CAD4TB (ver. 7) had the highest overall accuracy (73.8% specificity at 90% sensitivity), was significantly more specific than other algorithms, and achieved the minimum WHO target accuracy for a TB triage test. Additionally, the CAD4TB performance was stable regardless of diabetes status, suggesting its added value in TB-endemic regions with higher diabetes prevalence. (preprint) (Worodria, W. et al. 2024)

A study in Bangladesh showed CAD4TB (ver.7) achieved a sensitivity of 91.5% (fixed specificity at 70%) and a specificity of 72.8% (fixed sensitivity at 90%), concluding that CAD4TB performed better than human readers and meeting the WHO TPP values (Qin, ZZ. et al. PLOS DIGITAL HEALTH, 2022).

Another study in Kenya demonstrated that CAD4TB (ver. 6) achieved a sensitivity of 98.3% (fixed specificity at 70%) and a specificity of 90.4 % (fixed sensitivity at 90%) and met the optimal TPP (sensitivity of 95% and specificity of 80%) for TB community screening (Mungai, B. et al. PLOS GLOBAL PUBLIC HEALTH, 2022).

### Publication Highlights



CAD4TB met WHO's Target Product Profile values  
(Sensitivity > 0.90 and Specificity > 0.70)

Multiple studies



Performs on par or significantly better than human  
expert readers.

Multiple studies

<sup>1</sup> WHO consolidated guidelines on tuberculosis. Module 2: screening – systematic screening for tuberculosis disease. Geneva: World Health Organization; 2021. Licence: CC BY-NC-SA 3.0 IGO. <https://www.who.int/publications/i/item/9789240022676>



Daily screening yields 2.5 times higher with CAD4TB in comparison to screening without.

Murphy et al. Nature Scientific Reports. 2020.



Case detected 22% more with CAD4TB being used as a second CXR reader.

Philipsen et al. IJTLD. 2019.



With high throughput scenarios, the per-screen cost for CAD4TB is 73% lower than a radiologist for ACF and 61% lower for facility-based screening.

Bashir et al. PLOS One. 2022



Xpert tests used 75% savings when using CAD4TB at a sensitivity of 79.7%. More than 50% saving at a sensitivity CAD4TB above 90%.

Qin et al. The Lancet Digital Health, 2021.



TB & Diabetes screening: High accuracy with CAD4TB as a TB screening tool among people with diabetes.

Habib et al. Nature Scientific Reports. 2020.



TB & HIV: added-value with universal HIV screening increased timelines and completeness of HIV/TB diagnosis.

MacPherson et al. PLOS Medicine. 2021.



CAD has the potential to be a useful additional diagnostic tool for paediatric tuberculosis.

Palmer et al. PLOS Global Public Health. 2023.

A complete list of CAD publications is included in the submission.

## Performance

CAD4TB had the highest overall accuracy (73.8% specificity at 90% sensitivity), was significantly more specific than other algorithms, and achieved the minimum WHO target accuracy for a TB triage test.	[India, Madagascar, South Africa, Tanzania, the Philippines, South Africa, and Vietnam: <a href="#">Worodria W et al. 2024</a> ] <b>(preprint)</b>
The AUC (95% CI) of CAD4TB against the microbiological reference standard (Xpert Ultra and/or sputum culture positivity) was 0.90 (0.82-0.97).	[South Africa, Zambia, Zimbabwe: <a href="#">Scott A et al. International Journal of Infectious Diseases. 2024</a> ]
As threshold determination must be context-specific, our analytically-straightforward approach should be adopted to leverage prevalence surveys for CAD threshold determination in other settings with a comparable proportion of eligible but not tested participants.	[Lesotho: <a href="#">Vanobberghen F et al. ERJ Open Res. 2023</a> ]

The overall performance between CAD4TB v5, v6 and v7 was similar, but the distribution of scores across the 100-point scale varied greatly across the three versions. To accommodate intra-version variation, screening programmes need to select new triaging thresholds for each new software update.	[South Africa: <a href="#">Fehr J et al. IJTLD. 2023</a> ]
The newer version (CAD4TB 7) significantly outperformed the predecessor (CAD4TB 6), performing better than human readers and meeting WHO TPP values.	[Bangladesh: <a href="#">Qin ZZ et al. PLOS Digital Health. 2022</a> ]
CAD4TB v7 and two other CAD software solutions emerged from this evaluation as excellent alternatives for human CXR interpretation, performing on par with the Expert Reader and significantly better than the Intermediate Reader.	[Vietnam: <a href="#">Codlin A et al. Nature Scientific Report. 2021</a> ]
CAD4TB performs significantly better than human readers.	[Qin ZZ et al. <a href="#">The Lancet Digital Health. 2021</a> ] [Tanzania: <a href="#">Breuninger et al. PLOS One. 2014</a> ]
CAD4TB is an accurate tool for community-based TB screening for prevalence surveys in Kenya. CAD4TB 6 met the optimal WHO TPP.	[Kenya: <a href="#">Mungai B et al. PLOS Global Public Health. 2022</a> ]
CAD4TB is on par with human expert readers.	[Vietnam: <a href="#">Codlin A et al., Nature Scientific Report, 2021</a> ], [South Africa: <a href="#">Fehr et al., npj digital medicine, 2021</a> ], [G. Tavaziva et al., <a href="#">Clinical Infectious Diseases, 2021</a> ], [Pakistan: <a href="#">Murphy et al., Nature Scientific Reports, 2020</a> ] [Philippines: <a href="#">Philipsen et al., IJTLD, 2019</a> ], [London: <a href="#">Melendez et al., IJTLD, 2018</a> ], [Zambia: <a href="#">Melendez et al., IJTLD, 2017</a> ], [South Africa: <a href="#">Hogeweg et al., IEEE Trans Med Imaging, 2015</a> ], [Tanzania: <a href="#">Steiner et al. Public Health Action. 2015</a> ]
CAD4TB met the WHO's TPP at 90% sensitivity.	[Qin ZZ et al. <a href="#">The Lancet Digital Health. 2021</a> ]
In a community-based multi-disease screening survey in an HIV-endemic rural area, CAD4TB achieved comparable sensitivity and specificity to the radiologist. CAD4TB has the potential to replace radiologists for triaging CXRs in prevalence surveys.	[South Africa: <a href="#">Fehr et al. npj digital medicine. 2021</a> ]
In Tanzania, CAD4TB performance was significantly better than a clinical officer's.	[Tanzania: <a href="#">Breuninger et al. PLOS One. 2014</a> ]

## Efficiency & Effectiveness

<p>CAD has the potential to be a useful and cost-effective screening tool for TB in a resource-poor HIV-endemic African setting, assisting active case finding strategies to break the TB transmission cycle.</p>	<p>[Kenya, South Africa, Zambia: <a href="#">Scott AJ et al. Open Forum Infectious Diseases. 2024</a>]</p>
<p>In the context of community-based ACF in endemic TB/HIV settings, using POC Xpert and x-ray screening with CAD analysis is both feasible and had high diagnostic yield for TB and COVID-19.</p>	<p>[South Africa, Zambia, Zimbabwe: <a href="#">Scott AJ et al., International Journal of Infectious Diseases. 2024</a>]</p>
<p>Ensuring CXR for all can be operationally challenging in a programmatic setting, however, the combined use of digital X-ray and CAD4TB in this project replaced the human processing and interpretation of X-ray and contributed to the high CXR coverage.</p> <p>The proportion of people diagnosed with TB who had symptoms was very small. This suggested that most of the detected TB cases (&gt;80%) were asymptomatic and were captured by the contribution of CXR.</p>	<p>[Papua New Guinea: <a href="#">Dakulala et al., BMC Public Health. 2024</a>]</p>
<p>An experimental calibration method achieved a viable CAD threshold for testing. High CAD scores can identify subclinical TB and those at risk of progression to bacteriologically confirmed TB disease in the near term.</p> <p>To maximize feasibility, acceptability, quality, and yield of ACF over time, continuous review and re-calibration of outreach strategies, algorithm, and procedures were essential.</p>	<p>[Nigeria: <a href="#">Eneogu RA et al. PLOS Global Public Health. 2024</a>]</p>
<p>Delft Light portable digital X-ray and CAD4TB in parallel with the WHO 4-symptom screen achieved a lower pre-diagnostic loss of presumptive TB cases.</p> <p>The W4SS + portable X-ray with CAD screen-based intervention was more efficient for TB case yield: four times the number of TB cases yield than symptom-only-based screening intervention.</p>	<p>[Nigeria: <a href="#">Babayi et al. Public Health Action. 2023</a>]</p>
<p>CAD may provide viable options to increase TB detection, especially in low-resource areas where there may be no available expert radiologists.</p>	<p>[Gelaw SM et al. <a href="#">PLOS Global Public Health. 2023</a>]</p>
<p>Ultra-portable X-rays with CAD were overall well received to decentralise radiological assessment for TB.</p>	<p>[Qin ZZ et al. <a href="#">PLOS ONE. 2023</a>]</p>
<p>The per-screen costs for the two CAD software programs with a perpetual licensing costing structure are considerably lower than the cost with radiologists for high throughput in ACF scenarios.</p>	<p>[Pakistan: <a href="#">Bashir S et al. PLOS ONE. 2022</a>]</p>

With high throughput scenarios, the per-screen cost for CAD4TB is 73% lower than a radiologist for ACF and 61% lower for facility-based screening.	
TB screening using Delft Light Backpack X-ray and CAD4TB during community-based ACF in hard-to-reach Niger Delta communities of Nigeria showed a high TB prevalence among participants. Nationwide deployment of the instrument in hard-to-reach areas is recommended.	<a href="#">[Nigeria: Odume B et al. Public Health Action. 2022]</a>
CAD4TB v7 had a significantly higher AUC than v6, 0.903 (95% CI: 0.897–0.908) compared to 0.823 (0.816–0.830).  CAD4TB 7 results in a steep initial increase of Xpert test saved. Until a threshold of 75, greater numbers of Xpert tests can be saved using CAD4TB v7.	<a href="#">[Bangladesh: Qin ZZ et al. PLOS Digital Health. 2022]</a>
A combination of chest X-ray analysis by CAD4TB and symptomatology is of immense value in screening a large population at risk in a developing, high-burden country. It is significantly a more effective tool for screening and early diagnosis of TB in individuals who would otherwise go undiagnosed.	<a href="#">[Pakistan: Nishtar et al. Pak J Med Sci. 2022]</a>
Using CAD4TB can save molecular tests such as Xpert tests.	<a href="#">[Qin ZZ et al. The Lancet Digital Health. 2021]</a> , <a href="#">[Pakistan: Habib et al. Scientific Reports. 2020]</a> , <a href="#">[Pakistan: Zaidi et al. Nature Scientific Reports. 2018]</a> , <a href="#">[Bangladesh: Rahman T et al. European Respiratory Journal. 2017]</a>
All five AI algorithms (incl. CAD4TB) significantly outperformed the radiologists and reduced the number of Xpert tests required by 50% while maintaining a sensitivity above 90%.	<a href="#">[Qin ZZ et al. The Lancet Digital Health. 2021]</a>
Triage using AI-based CXR interpretation can be cost-effective and even cost-saving relative to standard practice among persons with suspected TB in Pakistan. These results are likely applicable to other low-income, high-TB-burden settings.	<a href="#">[Pakistan: Nsengiyumva et al. Open Forum Infectious Diseases. 2021]</a>
Community-based screening with innovative activities, comprising sensitive screening and diagnostic tools, including a digital X-ray and CAD4TB, effectively improves TB case detection.	<a href="#">[Pakistan; A. Wali et al. Public Health Action. 2021]</a>
CAD4TB, as a triage tool to pre-screen people with diabetes for Xpert testing, can improve case-detection while potentially reducing program costs with more rational use of expensive molecular tests.	<a href="#">[Pakistan; Habib et al. Nature Scientific Reports. 2020]</a>

In Pakistan, the cost per screened subject using CAD4TB v6 is almost half in comparison to screening without, while almost 2.5 times higher daily throughput.	<a href="#">[Pakistan; Murphy et al. Nature Scientific Reports. 2020]</a>
Used as a pre-screening system in TB-endemic regions, CAD4TB allows for the testing of much larger numbers of subjects at a fraction of the cost.	<a href="#">[Pakistan; Murphy et al. Nature Scientific Reports. 2020]</a>
In the Philippines, CAD4TB as a second reader showed a 22% increase in TB detection and could help detect additional TB-positive cases that were initially missed at limited additional costs.	<a href="#">[Philippines: Philipsen et al. IJTL. 2019]</a>
In Pakistan, CAD4TB, as a triage tool, could minimise the use of expensive molecular tests in low-resource settings.	<a href="#">[Pakistan; Zaidi et al. Nature Scientific Reports. 2018]</a>
Reduces cost per identified TB case and the cost per screened subject by almost half.	<a href="#">[Philippines: Philipsen et al., Nature Scientific Reports. 2015]</a>
Useful where trained human readers are scarce.	<a href="#">[Zambia: Muyoyeta et al. PLOS One. 2014]</a>

## Paediatric TB

The performance of CAD4TB v7 to identify TB in children (<13 years) significantly improved after fine-tuning it with a set of well-characterised paediatric chest x-rays. CAD has the potential to be a useful additional diagnostic tool for paediatric tuberculosis.	<a href="#">[South Africa: Palmer et al. PLOS Glob Public Health. 2023]</a>
--	---

## Key and Vulnerable Populations

### TB / HIV

Urine SILVAMP-LAM/Determine-LAM plus dCXR-CAD4TB diagnostics identified more hospitalised PHIV with TB than usual care.	<a href="#">[Malawi: Burke, RM et al. Clinical Infectious Diseases. 2024]</a>
In a community-based multi-disease screening survey in an HIV-endemic rural area, CAD4TB achieved comparable sensitivity and specificity to the radiologist. HIV serostatus did not impact CAD4TB's performance.	<a href="#">[South Africa: Fehr et al. npj digital medicine. 2021]</a>
Digital CXR using CAD4TB with universal HIV screening significantly increased timelines and completeness of HIV and TB diagnosis.	<a href="#">[Malawi: MacPherson et al. PLOS MEDICINE. 2021]</a>

Combined use of mobile X-ray and CAD4TB to identify TB among high-risk groups to improve early TB diagnosis can deliver a major contribution to TB control in Romania.	<a href="#">[Romania: Mahler B et al. BMJ Open. 2021]</a>
In an HIV-negative population, CAD4TB v6 met WHO-recommended minimal accuracy for pulmonary tuberculosis triage tests.	<a href="#">[Pakistan: Khan FA et al. The Lancet Digital Health. 2020]</a>

## Diabetes

CAD4TB performance was stable regardless of diabetes status.	[India, Madagascar, South Africa, Tanzania, the Philippines, South Africa, and Vietnam: <a href="#">Worodria W et al. 2024</a> ]
CAD4TB offers good diagnostic accuracy as triage for TB screening among diabetes patients.	<a href="#">[Pakistan: Habib et al. Nature Scientific Reports. 2020]</a>
CAD4TB with X-ray systems allowed not only rapid and systematic triage to Xpert testing but also found quantitatively more TB-like abnormalities in those with Diabetes Mellitus.	<a href="#">[Bangladesh: Paul KK et al., Science Direct. 2020]</a>
In Indonesia, CAD4TB has the potential as a triage tool for TB screening in people living with diabetes, thereby significantly reducing the need for microbiological examination.	<a href="#">[Indonesia: Koesoemadinata RC et al. IJTLD. 2018]</a>

## Screening in Prison

Screening by mobile X-ray systems with automated interpretation could reduce the number of confirmatory tests required and enable screening to be more rapid in high-burden TB settings while still maintaining sufficient sensitivity.	[Brazil: <a href="#">Soares et al., The Lancet Regional Health – Americas. 2023</a> ]
The inclusion of digital CXR in systematic TB screening detected additional TB cases among inmates that would otherwise have been missed, and using CAD4TB may also improve the performance of the screening algorithm.	[South Africa: <a href="#">Kim et al., IJTLD. 2020</a> ]
High uptake of new screening tools, such as digital X-rays with CAD4TB, may be particularly feasible, reliable, and highly acceptable in prison settings.	[Pakistan: <a href="#">Wali et al. BMC Public Health. 2019</a> ]
CAD4TB reliably evaluates CXRs from a mostly asymptomatic prison population, with a performance comparable to local readers in Tanzania.	[Tanzania: <a href="#">Steiner et al., Public Health Action. 2015</a> ]

## Non-TB Abnormalities

CAD4TB has the potential to simultaneously provide information on other non-TB abnormalities that might be of clinical relevance in communities alongside TB.

CAD can be useful to LMICs where there is no routine screening for non-TB abnormalities, and there is often a shortage of qualified radiologists.

[Zambia, South Africa: [Ngosa D. et al. BMC Infect Dis. 2023](#)]

## Silicosis

CAD4Silicosis achieved a high AUC of 0,926 and 0,903 against readers 1 and 2 (silicosis classified as ILO $\geq$ 1/1). On the 90% sensitivity criterion, a specificity of 82.6% against reader 1 and 74.9% against reader 2.

[South Africa: [Ehrlich R et al. Int. J. Environ. Res. Public Health. 2022](#)]